


**Supplementary Figure 1. The transfection efficiency of CCDC65 plasmid, lentivirus and siRNA.** (A-B) The expression of CCDC65 in GC cells treated with CCDC65 lentivirus and plasmid. (C-D) The fluorescence of GFP was examined in LV-CCDC65 and LV-NC GC cells (scale bar: 50  $\mu$ m). (E) The expression of CCDC65 in GC cells treated with CCDC65 siRNA and control was examined by QPCR.




**Supplementary Figure 2. The 130-484aa domain of CCDC65 suppressed cell proliferation and metastasis in GC. (A)**


MTT assays of two CCDC65 truncated mutants and control group in GC cells. All data are presented as the mean  $\pm$  SD. **(B)**

Transwell and boyden assay evaluated the migration and invasion of gastric cancer cells treated with two CCDC65 truncated mutants and control vector. Student's t test. **(C)** Expression levels of ENO1 and p-AKT1 (ser473) were detected following


transfection with CCDC65 truncated mutants and control vector. GAPDH was used as a loading control.




**Supplementary Figure 3. CCDC65 recruits FBXW7 to regulate ENO1 protein stability.** (A-B) WB was used to detect the effects of DMSO or MG132 treatment and CHX treatment for different duration on the stability of ENO1 protein in the CCDC65 and CCDC65+si-FBXW7 groups. (C) GST-CCDC65 interacts with His-FBXW7 in vitro by GST pull-down assay. (D) Co-IP was conducted to identify the function of CCDC65 knockdown on the interaction among ENO1, FBXW7 and ubiquitin in GC cells.



**Supplementary Figure 4. Ectopic expression of ENO1 mitigates CCDC65 suppression of GC proliferation and metastasis. (A)** MTT assays were used to detect the viability of GC cells in control, control+ENO1, CCDC65+control, CCDC65+ENO1 groups. **(B)** Transwell and Boyden assay were carried out in control, control+ENO1, CCDC65+ control, CCDC65+ENO1 groups. **(C)** The expression of ENO1, p-AKT1 (Ser473), N-cadherin, E-cadherin, Vimentin, CCND1 and P21 in control, control+ENO1, CCDC65+ control, CCDC65+ENO1 groups.



**Supplementary Figure 5. Metformin potentiates CCDC65-mediated suppression of proliferation, metastasis and ENO1-AKT1 signal in GC. (A-B)** MTT, transwell and boyden assay were carried out in si-NC, si-CCDC65, metformin+ si-NC, metformin+ si-CCDC65 groups. **(C)** The weight of mice treated with normal saline (NS), metformin (met) alone or metformin in combination with si-CCDC65 was collected every 3 days.



**Supplementary Figure 6. The percent of ki67 staining of xenograft. (A)** The percent of ki67 staining in xenograft tumors originating from LV-NC and LV-CCDC65 cells. **(B)** The percent of ki67 staining in xenograft tumors originating from si-NC and si-CCDC65 cells.

**Supplementary Table 1: Association of CCDC65 expression with gastric cancer clinicopathological characteristics.**

| Parameters            | Category | NO. | CCDC65 expression |           | P      |
|-----------------------|----------|-----|-------------------|-----------|--------|
|                       |          |     | Low (97)          | High (90) |        |
| Age (years)           | < 60     | 64  | 32                | 32        | 0.7118 |
|                       | ≥ 60     | 123 | 65                | 58        |        |
| Gender                | Male     | 133 | 68                | 65        | 0.7494 |
|                       | Female   | 54  | 29                | 25        |        |
| T stage               | T1+T2    | 29  | 7                 | 22        | 0.0012 |
|                       | T3+T4    | 158 | 90                | 68        |        |
| N stage               | N0+N1    | 93  | 35                | 58        | 0.0001 |
|                       | N2+N3    | 84  | 62                | 22        |        |
| M stage               | M0       | 157 | 75                | 82        | 0.0157 |
|                       | M1       | 30  | 22                | 8         |        |
| Clinical stage (AJCC) | I+II     | 64  | 22                | 42        | 0.0006 |
|                       | III+IV   | 123 | 75                | 48        |        |

**Supplementary Table 2: The sequences used in this study.**

| gene   | No. | target sequences    |
|--------|-----|---------------------|
| CCDC65 | 1   | GACTCACCTGGAAAGTAA  |
|        | 2   | CCAGAGAACTTCATAAGGA |
|        | 3   | CAAGGAGTTTGAGACAGAA |
| ENO1   | 1   | AGTCCTTCATCAAGGACTA |
|        | 2   | GCTGCTGAAGACTGCTATT |
|        | 3   | GGAAGTATGACCTGGACTT |
| FBXW7  | 1   | GGAACCCAAAGACCTGCTA |
|        | 2   | GTTAGTGGTTCTGATGACA |
|        | 3   | GCGTTGTATGCATCTTCAT |

**Supplementary Table 3: The primers used in this study.**

| Primers name |         | Sequence (5'-3')           |
|--------------|---------|----------------------------|
| CCDC65       | Forward | AGAAGTGCCTTCGGGAAGTC       |
|              | Reverse | CTCGTTCAAATGTTTGGCTGAG     |
| ENO1         | Forward | GCCGTGAACGAGAAGTCCTG       |
|              | Reverse | ACGCCTGAAGAGACTCGGT        |
| GAPDH        | Forward | TGCACCACCAACTGCTTA         |
|              | Reverse | GGATGCAGGGATGATGTTC        |
| ACTB         | Forward | GACCTGACTGACTACCTCATGAAGAT |
|              | Reverse | GTCACACTTCATGATGGAGTTGAAGG |

**Supplementary Table 4. A list of antibodies used for Western blot, IHC, IF, Co-IP.**

| antibody              | Company     | Cat.No     | Mol weight (kDa) | Dilution                                 |
|-----------------------|-------------|------------|------------------|--|
| CCDC65                | Proteintech | 24376-1-AP | 57               | 1:1000 (WB)<br>5 µg (IP)                 |
| CCDC65                | Abcam       | ab122482   | 57               | 1:50 (IHC)<br>1:100 (IF)                 |
| ENO1                  | Proteintech | 11204-1-AP | 47               | 1:5000 (WB)<br>1:250 (IHC)<br>1:100 (IF) |
| ENO1                  | CST         | 3810       | 47               | 1:50 (IP)                                |
| AKT1                  | SAB         | 48884      | 60               | 1:1000 (WB)<br>1:100 (IF)                |
| Phospho-AKT1 (Ser473) | CST         | 9018       | 60               | 1:3000 (WB)<br>1:400 (IHC)               |
| Nca                   | Proteintech | 66219-1-Ig | 130              | 1:2000 (WB)                              |
| Eca                   | Proteintech | 60335-1-Ig | 120              | 1:2000 (WB)                              |
| vimentin              | Proteintech | 10366-1-AP | 54               | 1:1000 (WB)                              |
| CCND1                 | Proteintech | 60186-1-Ig | 36               | 1:5000 (WB)                              |
| P21                   | CST         | 2947       | 21               | 1:1000 (WB)                              |
| FBXW7                 | Proteintech | 28424-1-AP | 100-110          | 1:1000 (WB)<br>5 µg (IP)                 |
| ubiquitin             | Proteintech | 10201      | 8, 24            | 1:500 (WB)                               |
| PCNA                  | Proteintech | 10205-2-AP | 36-38            | 1:500 (IHC)                              |
| Ki67                  | Proteintech | 27309-1-AP | -                | 1:8000 (IHC)                             |
| FLAG                  | Sigma       | F1804      | -                | 1:1000 (WB)<br>5 µg (IP)                 |
| MYC-tag               | Proteintech | 16286-1-AP | -                | 1:1000 (WB)<br>5 µg (IP)                 |
| HA-tag                | Proteintech | 51064-2-AP | -                | 1:3000 (WB)<br>5µg (IP)                  |
| GST                   | Proteintech | 10000-0-AP | -                | 1:1000 (WB)<br>5 µg (IP)                 |
| His                   | Proteintech | 66005-1-Ig | -                | 1:1000 (WB)<br>5 µg (IP)                 |
| GAPDH                 | Proteintech | 60004-1-Ig | 36               | 1:5000 (WB)                              |
| ACTB                  | Proteintech | 66009-1-Ig | 42               | 1:5000 (WB)                              |

WB: western blot

IHC: immunohistochemistry

IF: immunofluorescence

Co-IP: co-immunoprecipitation